



PORK DIPSTICK TEST

Doubt?

#porkscreening #DIY #fast #let's #self-test #haltest

Fast & on-the-spot screening of PORK TRACES

Accurate . Fast . Portable .



Raw Meat



Processed Meat / Canned Food



Gelatin



Surface Swab

Integrated Solution for Halal Food Testing

www.7foodpillars.com

Find us on: facebook.


facebook.com/7FoodPillars

PORK Free?



#fastresult #10-20minutes

Direction Of Use



Fast

Result within
10-20 minutes



Visible Result

If two red lines
appear at the strip,
food sample is
contaminated with
pork



Accurate



On-the-spot

Immediate
screening
anytime,
anywhere



Easy to Use

Simple protocol, can
be operated by
non-technical users



Travel Friendly

1 test strip per pack,
light and easy to
keep and carry



Validated by ISO 17025 Accredited Laboratory

Performance of
the test strips
were validated by
ISO 17025 accredited
laboratory

Ordering Information

Product Code	Description	Pack Size
RHT01-01-001	HALTEST Pork Dipstick Test - Raw Meat	1 test/ kit
RHT01-02-001	HALTEST Pork Dipstick Test - Processed Meat/ Canned Food	1 test/ kit
RHT01-03-001	HALTEST Pork Dipstick Test - Gelatin	1 test/ kit
RHT01-04-001	HALTEST Pork Dipstick Test - Surface Swab	1 test/ kit